AS 200E

Instructions, Gebrauchsanweisung,
Instrucciones, Instructions, Istruzioni



Contents



Labels and Symbols	4
EN/English.....	7
Notice	7
Description	8
Setting Up the Scanner	15
Getting Started	20
Safety Guide	24
Troubleshooting.....	32
Regulatory Information.....	34
Technical Specifications	41
DE/Deutsch.....	43
Hinweis	43
Beschreibung	44
Einrichtung.....	51
Erste Schritte	56
Sicherheitsleitfaden	60
Fehlersuche	69
Regulatorische Informationen.....	71
Technische Daten	78
ES/ Español	80
Aviso	80
Descripción.....	81
Configuración del Escáner	88
Primeros Pasos.....	93

Guía de Seguridad	97
Solución de Problemas	105
Información Reglamentaria	107
Especificaciones Técnicas	114
FR/Français	116
Avis.....	116
Description	117
Configuration du Scanner	124
Pour Commencer	129
Guide de Sécurité	133
Dépannage	142
Informations Réglementaires	144
Spécifications Techniques	151
IT/Italiano	153
Avviso.....	153
Descrizione	154
Impostazione del Scanner	161
Come Iniziare.....	166
Guida alla Sicurezza	170
Risoluzione dei Problemi	179
Informazioni Normative.....	181
Specifiche Tecniche	188
Contact Information.....	190

Labels and Symbols

	EN Type BF applied part DE Anwendungsteil vom Typ BF ES Componente aplicado de tipo BF FR Partie appliquée de type BF IT Parte applicata di tipo BF
	EN Separate collection for waste electric and electronic equipment (WEEE) is required. DE Die getrennte Sammlung von Elektro- und Elektronikgeräten (WEEE) ist vorgeschrieben. ES Es necesario recoger por separado los residuos de aparatos eléctricos y electrónicos (RAEE). FR La collecte sélective des déchets d'équipements électriques et électroniques (DEEE) est obligatoire. IT È richiesta la raccolta differenziata dei rifiuti di apparecchiature elettriche ed elettroniche (RAEE).
	EN Manufacturer DE Hersteller ES Fabricante FR Fabricant IT Fabbrikante
	EN Date of Manufacture DE Herstellungsdatum ES Fecha de manufactura FR Date de fabrication IT Data di Produzione
	EN Refer to instruction manual/booklet DE Siehe Gebrauchsinformation ES Consultar el manual de instrucciones/manual FR Se reporter au manuel/notice d'utilisation IT Consultare le istruzioni del Manuale d'Uso
	EN Direct current DE Gleichstrom ES Corriente continua FR Courant continu IT Corrente continua
	EN Sterilizable in a steam sterilizer (autoclave) at 134°C DE Sterilisierbar in einem Dampfsterilisator (Autoklav) bei 134°C ES Esterilizable en esterilizador de vapor (autoclave) a 134 °C FR Stérilisable dans un stérilisateur à vapeur (autoclave) à 134°C IT Sterilizzabile in uno sterilizzatore a vapore (autoclave) a 134°C

	EN Model Number DE Modellnummer ES Número de modelo FR Numéro de modèle IT Numero di modello
	EN Serial number DE Seriennummer ES Número de serie FR Numéro de série IT Numero di serie
	EN Rechargeable battery DE Wiederaufladbare Batterie ES Batería recargable FR Batterie rechargeable IT Batteria ricaricabile
	EN CE marking applicable for European Union DE CE-Kennzeichnung gilt für die Europäische Union ES Marca CE aplicable para la Unión Europea FR Marquage CE applicable à l'Union européenne IT Marchio CE applicabile all'Unione Europea
	EN Authorized Representative in the European Community DE Bevollmächtigter in der Europäischen Gemeinschaft ES Representante autorizado en la Comunidad Europea FR Mandataire Européen IT Rappresentante autorizzato per la Comunità Europea
	EN Authorized Representative in Switzerland DE Autorisierter Vertreter in der Schweiz ES Representante autorizado en Suiza FR Représentant autorisé en Suisse IT Rappresentante autorizzato in Svizzera
	EN Medical Device DE Medizinisches Gerät ES Productos sanitarios FR Dispositif médical IT Dispositivo medico
	EN Atmospheric pressure limitation DE Atmosphärische Druckbegrenzung ES Limitación de la presión atmosférica FR Limitation de la pression atmosphérique IT Limitazione della pressione atmosferica
	EN Fragile, handle with care DE Zerbrechlich, mit Vorsicht zu behandeln ES Frágil, manipular con cuidado

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EN/English

Notice

This document includes information on the safety instructions, regulatory information, Technical Specifications and installation and operation instructions of the device. We recommend that you thoroughly familiarize yourself with this document to make the most effective use of your system.

The content of this document may change over time. Alliedstar and its subsidiaries are not responsible for any errors in this document or any incidental damages related to its provision, performance, or use. Reproduction of any part of this publication is prohibited without Alliedstar's permission.

Conventions

The following special notices highlight important information or potential risks to personnel or equipment.

	WARNING: Indicates safety instructions that must be followed precisely to prevent injury to yourself or others.
	CAUTION: Warns of conditions that could cause severe damage or operational issues.
	NOTE: Provides additional information and helpful tips.

WARNING: We recommend that you consult this document before using the device.

All trademarks and registered trademarks are the property of their respective holders.

The device is intended for professional use only.

U.S. Federal law restricts this device to sale by or on the order of a dentist.

If any serious incident occurs in relation to the device, the user must report it to Alliedstar and to the competent authority of its Member State in the European Union.

Description

Intended Use

The **AS 200E Intraoral Scanner** (herein after referred as ‘**scanner**’) is a digital optical scanning device used to obtain digital impressions of hard and soft tissues such as teeth, gingival, and mucous membrane by means of oral scanning, for oral restoration and orthodontic treatment of malocclusion.

The scanner could be used for both adults and children in clinical practice.

Contraindications

None


Clinical Benefits and Performance Characteristics

The scanner benefits dental practices by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

The scanner is designed to acquire 3D models in the followings:

- Upper jaw
- Lower jaw
- Buccal bite registration

Nomenclature

[1] Handpiece
[2] Scanner Tip
[3] Scan Button ● Press once to start scanning, press again to stop scanning
[4] Mode Indicator
 Upper jaw scan mode
 Lower jaw scan mode
 Buccal bite registration mode
[5] Mode Switch Button ● Press this button to switch between different modes
[6] Status Indicator ● Slowly flashing: starting up ● Rapidly flashing: preparing network ● Breathing: ready to be connected to the software ● Solid on: Connected to the software and selected as active scanner
[7] Battery Indicator In Using: ● Solid green: battery level higher than 50% ● Solid yellow: battery level between 20% and 50% ● Flashing yellow: battery level lower than 20%

In Charging: ● Breathing green: battery level higher than 50% ● Breathing yellow: battery level lower than 50% ● Solid green: fully charged
[8] Lens Window
[9] Charging Holder
[10] USB Cable (1m length)

Computer System Requirements


For the computer system requirements, refer to the 'Technical Specifications' section


WARNING: It is MANDATORY to check that your computer system configuration is compatible with the computer system requirements for the accompanying software.

Charging the Battery

A fully charged battery can support up to 2 hours of scan time. It is recommended to keep the scanner in the powered charging holder when it is not in use. When the Battery Indicator flashes yellow, it means the battery is low and needs to be charged in time. To charge, the scanner needs to be inserted into the Charging Holder.


The Charging Holder must be connected to a USB Type C power source to be functional.


You can also continue to use the scanner by replacing the battery, follow these steps to replace the battery.

- Turn the scanner battery cover with the battery cover turner in the wall mounting kit or a coin to unlock.


- Take out the low battery and insert a charged one.


WARNING: Please use the battery specified by Alliedstar (Li-18650-3.6V 3500mAh-PCM-Cap)

- Close the scanner battery cover by turning it to lock.


- The charging holder holds a spare battery. Open the spare battery cover by turning it to unlock and you will be able to take the spare battery out.
- **NOTE:** The spare batteries inside the charging holder are not charged when the scanner is in the charging holder (the spare batteries inside the charging holder are only charged when the scanner is removed from the charging holder).

NOTE: Please fully charge the battery before use.


Wireless Connection

The scanner connects to the software through a wireless connection. Before using the scanner, a compatible wireless adapter needs to be connected to the computer that runs the software.

NOTE: Expand the wireless adapter to maximize the network signal strength. The wireless card product model number can be viewed at the bottom of <https://www.allied-star.com/product/as-200e>.


After the compatible wireless adapter is connected, select the scanner you want to connect from the Scanner List window at the bottom right of the software and click the connect button.

	
Indicates a wireless scanner	Indicates a wired scanner

Charging Holder Overview

Place the scanner in the Charging Holder when you are not using it.


NOTE: The scanner will automatically shut down if it is idle for more than 3 minutes. If you need to use it again, press any button on the scanner, or take it out of the Charging Holder.

Setting Up the Scanner

Setting Up the Scanner

To set up the scanner, follow these steps:

- Visit <https://support.allied-star.com/download> and download the installation file and instructions according to the product model.
- Double-click the software installation file.
- Choose a language from the drop-down list and click Ok to install.
- Follow the instructions on the screen to complete the installation.
- Remove the lens window protective cover on the Lens Window.
- Firmly slide the Scanner Tip onto the end of the scanner.

- Short press any button to start the scanner.
- After the scanner is powered on, open the software and connect the scanner through the wireless network.


Using the Scanner Charging Holder

The Charging Holder can be used standing on a desktop or mounted on the wall.


Installing the Desktop Charging Holder

To place the Charging Holder on a desktop, follow these steps:


- Select a clean surface area.
- Insert one end of the USB cable into the bottom of the charging holder.


- Insert the charger with the connected USB cable into the stand.

- Insert the USB cable into the cable clips.

- Pass the other end of the USB cable through the hole reserved on the stand.


- Insert the scanner into the desktop charging holder.

Installing the Wall Mount Charging Holder

To mount the Charging Holder on the wall, follow these steps:


- Select an area that you can access easily.


- Insert the screws through the holes in the holder rail and wall mount pad in turn to affix it to a solid surface.
- Insert one end of the USB cable into the bottom of the charging holder.

- Insert the Charging Holder into the holder rail.

- Insert the scanner into the Wall Mount Charging Holder.

CAUTION: If the Wall Mount Charging Holder is not properly installed, there is a risk that the Charging Holder can fall off the wall, resulting in damage to the scanner.

Getting Started

Accessing the Software User Interface

To access the software user interface, follow these steps:

- Double-click the software icon on the desktop.
- The software will automatically pop up the login window.
- If you don't have an account, click **Create an account** to register the organization, and complete the email verification.
- Type your account information in the username and password fields, click the **Log in** button.
- The **Cases** page will be shown up.
- Click the **New Case** button, then input patient name and click the **OK** button.
- Click the **Scan** button.
- If the scanner is not activated, the device activation dialog will be displayed. Follow the instructions on the screen to complete device activation.
- Click the Option menu button  and select **Preferences**.
- Customize the configuration options.
- You can now start acquiring 3D models.

Audio Prompts

Scanning Sound

When the **Enable scanning sound** option is enabled, your computer will play a continuous sound when you are successfully scanning. If the sound stops, it means the scan has stopped. If you need to continue, please return to the previous scanning area until the scanner resumes scanning and your computer plays a continuous sound. When a buccal bite registration is successfully scanned, your computer will also play a short sound.

Warning Sound

When the **Enable warning sound** option is enabled, if the cumulative scanning time of the current case exceeds the recommended threshold and your computer may not be able to maintain the peak scanning performance, your computer will play a short warning sound.

NOTE: If your computer does not have speakers, these options will not take effect.

Preparing the Teeth

- If there is a preparation area, retract the gingiva by gingival retraction cords. And extract the cords just before scanning the preparation.
- Before starting the scan, dry the teeth thoroughly.
- During the scan, re-dry the teeth moderately.

Preparing the Scanner


The Scanner Tip attached to the Handpiece provides sanitary shield for patients. Always disinfect the body of the scanner and perform sterilization on the Scanner Tip after each use.

CAUTION: Scanner Tips are shipped non-sterile. You must sterilize them before the first use.


CAUTION: Avoid any liquid from leaking into the air outlet near the scanner lens window or the air inlet at the rear of the Handpiece, otherwise the Handpiece may be damaged.

To prepare the scanner, follow these steps:

- Make sure the Lens Window at the base of the scanner is clean by wiping it with a moist, lint-free cloth or lens tissue.
- Slide the Scanner Tip onto the scanner as shown below.


Starting Scanning

To start scanning, place the Scanner Tip on the surface of the tooth to stabilize the scanner and press the Scan button. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-5mm from the teeth.

Scanning Approach

The recommended scanning method is to start with a molar, since it has greater details for easier identification. Change the scanning angle to less than 60 degrees during scanning to allow the surfaces to overlap, if the overlap is too small, the alignment may be lost.


Scanning Protocol

The recommended scanning protocol consists of 3 sweeps: occlusal, lingual and buccal to ensure good data coverage of all surfaces.

It is recommended to start the first sweep from the occlusal surface, and you should start with the first molar.

The second sweep can scan both the lingual and buccal sides, and the third scan covers the opposite side of the second sweep.

Safety Guide

Warnings and Safety Instructions

DANGER OF ELECTRIC SHOCK

This is an electrical unit. Do NOT expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.

CAUTION: All known residual risks, contraindications, or undesirable side effects are listed in this document. If any serious incident occurs in relation to the device, you must report it to Alliedstar and to the competent authority of your Member State in the European Union.

WARNING

Scanner

- Before using the scanner, you must read and understand all safety information.
- Use the scanner only within hospitals and professional healthcare facilities. Do not operate it near high-frequency surgical equipment or inside the RF-shielded room of an ME system for magnetic resonance imaging, where electromagnetic disturbances are intense.
- Before operating the scanner, inspect the scanner and its accessories for any rough surfaces, sharp edges, or protrusions that could pose a safety hazard.
- You are responsible for the scanner's operation and maintenance. Training is required to use the scanner safely.
- Do not place objects within the scanner's operational field.
- Ensure the scanner is turned off when not in use.
- Do not use the scanner in oxygen-rich environments. It is not intended for use with flammable anesthetics or agents.
- Avoid pulling or twisting the cable.
- Do not drop the scanner or its accessories.
- Do not attempt to sterilize the scanner.
- Keep the scanner away from water sprays, and do not immerse it in water or disinfectants.
- Avoid subjecting the scanner to heavy vibrations.

- Avoid exposing the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- Avoid staring directly at the LED emission window.
- When the Scanner Tip is removed, install the lens window protective cover to protect the Lens Window.
- Avoid removing the cover from any part of the scanner. The device has no user-serviceable components. If repairs are needed, reach out to a qualified Alliedstar service technician.
- Do not use any cables other than those provided with the scanner. Using different cables may damage the scanner and compromise its safety features and EMC performance.
- Keep equipment that does not comply with IEC-60601 standards at least 1.5 meters from the patient.
- If the equipment malfunctions, switch off, post an “Out of Service” notice, and contact a certified Alliedstar technician.
- The use of components, accessories, cables, or spare parts not specified or supplied by the manufacturer of this equipment can compromise the scanner’s safety features. This may lead to higher electromagnetic emissions, reduced electromagnetic immunity, and improper equipment functioning.
- Do not modify this equipment in any way.
- Do not connect additional power strips or extension cords to the system.
- The applied part may reach temperatures up to 43 °C (109.4 °F). To prevent overheating, avoid prolonged use.
- To isolate the Charging holder from power supply, unplug the USB connector from the USB port.
- DO NOT maintain or service this equipment while it is in use with the patient.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.
- Patients with oral mucosal disease, mental illness, severe respiratory disease, asthma, Parkinson's disease, hyperactivity disease are forbidden.
- Patients with moderate or severe opening limitation should use it with caution.

Computer

- Do NOT place any equipment which does not comply with IEC 60601-1 in the immediate vicinity of the patient. Leave at least 1.5 meters distance between the patient and the equipment.
- The scanner is only intended to be connected to a computer that is at least IEC 60950 / IEC 62368, or equivalent standards certified. Connecting the scanner to other equipment may be hazardous.
- Refer to your computer's installation guide for details on the data processing system, computer, and monitor. Ensure there is enough space around the computer to allow for proper ventilation.

- Place the screen so that light reflections from both internal and external sources are minimized to achieve the best image quality and ensure comfortable viewing.

Battery

- Do not dismantle, open, or shred secondary cells or batteries.
- Do not expose cells or batteries to heat or fire. Avoid storage in direct sunlight.
- Do not short-circuit a cell or a battery. Do not store cells or batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- Do not remove a cell or battery from its original packaging until required for use.
- Do not subject cells or batteries to mechanical shock.
- In the event of a cell leaking, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- Do not use any charger other than that specifically provided for use with the equipment. Refer to the manufacturer's instructions for proper charging instructions.
- Do not use any cell or battery which is not designed for use with the equipment.
- Always purchase the battery recommended by the device manufacturer for the equipment.
- Keep cells and batteries clean and dry.
- Wipe the cell or battery terminals with a clean, dry cloth if they become dirty.
- Do not leave a battery on prolonged charge when not in use.
- After extended periods of storage, it may be necessary to charge and discharge the cells or batteries several times to obtain maximum performance.
- Retain the original product literature for future reference.
- Use only the cell or battery in the application for which it was intended.
- When possible, remove the battery from the equipment when not in use.
- Please remove the battery using specified tools. Removing the battery by hand is not allowed.
- Dispose of the battery properly.

Disposal

This device contains materials and chemicals commonly used in electrical and electronic products. Improper disposal at the end of its life can lead to environmental pollution. Therefore, do not discard it with regular household waste. Please take it to a designated electronic waste disposal or recycling facility. For more information on how to dispose of electronic waste, please get in touch with your local authorities.

Dispose of the Scanner Tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional Scanner Tips, contact your dealer.

Cleaning, Disinfecting, Sterilizing

Perform the following maintenance activities on your scanner and accessories regularly. Regularly perform the maintenance tasks outlined below for your scanner and accessories.

Carefully follow the instructions when preparing the scanner for use to ensure the highest level of hygienic safety for patients.

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories.

After each patient:

- Clean and disinfect the scanner. See “[Clean and Disinfect the Scanner](#)”).
- Clean the Scanner Tip, and then perform Autoclave Sterilization (See “[Clean the Scanner Tip](#)” and “[Sterilize the Scanner Tip](#)”).

There are 3 product models of Scanner Tip:

Product Model	Size	UDI-DI	Manual Cleaning	134°C Sterilization
TP202	Large	(01)06973993441034	Yes	Yes
TP203	Small	(01)06973993441041	Yes	Yes
TP204	Large	(01)06973993441386	Yes	Yes

Clean and Disinfect the Scanner

General Warnings

WARNING

- When handling the scanner, read and follow all warnings and personal protection guidelines in the sanitizer's Safety Data Sheet (SDS).
- Gloves must be worn when cleaning and disinfecting the scanner.
- After each patient, disinfect the scanner using a recommended medium-level disinfectant wipe with tuberculocidal properties.
- Do not use disinfectants containing phenol or iodine-based substances, as they may damage the scanner's

surface coating.

- Do not place the scanner in a pressure steam sterilizer or immerse it in water or disinfectant solutions.
- Excessive liquid can harm the scanner.
- Avoid using cotton balls, cloth pieces, or paper towels soaked in disinfectant to sanitize the scanner.

Clean the Scanner

- If visible blood or body fluid is on the scanner, clean it before proceeding with disinfection.
- To clean the scanner, follow the below steps:
 - ◆ Dampen a lint-free cloth with warm water (do not soak).
 - ◆ Use the dampened lint-free cloth to wipe away blood or body fluids.

Disinfect the Scanner

The scanner must be thoroughly sanitized after use for each patient.

To properly disinfect the scanner, make sure to follow the disinfectant manufacturer's guidelines regarding the required contact time.

CAUTION: The scanner must be cleaned before disinfection if it is visibly soiled. Refer to "[Cleaning the Scanner](#)."

To disinfect the scanner, proceed with the following steps:

- Detach the Scanner Tip.
- Clear any visible debris (refer to "[Cleaning the Scanner](#)").
- Use an intermediate-level disinfectant wipe that is commercially available. Be sure to follow the manufacturer's instructions for the necessary contact time.

Recommended disinfectant wipes: CaviWipes

WARNING: Using an unapproved disinfectant may damage the scanner.

- Thoroughly wipe down all scanner surfaces, ensuring no liquid enters any gaps, air outlets, or pinholes.

WARNING: Do not rinse the scanner.

- Let it air dry thoroughly.

- Once dry, use a clean, lint-free cloth dampened with water to remove any remaining disinfectant from the scanner's surface.

Clean and Sterilize the Scanner Tip

WARNING

- Always wear gloves when handling a contaminated Scanner Tip.
- Review and follow the safety and personal protection instructions in the detergent manufacturer's Safety Data Sheet (SDS).
- Avoid soaking the Scanner Tip in disinfectant for extended periods.
- Ensure the scanner tip is completely dry before mounting it to the scanner.
- Do not use an ultrasonic cleaner to clean the Scanner Tip.
- Do not soak the Scanner Tip with alcohol-based disinfectants.

Clean the Scanner Tip

To manually cleaning the Scanner Tip, follow these steps:

- Rinse the scanner tip under running water to remove excess debris (approximately 2 minutes).
- Gently brush all surfaces with an enzymatic detergent solution, like Metrex EmPower.
- Rinse thoroughly under fresh, running water for a full 2 minutes.
- Inspect the scanner tip. If it remains dirty, go through the steps above again.
- Clean the mirror inside the Scanner Tip using a lint-free cloth or lens tissue to remove dust.
- Dry the scanner tip's mirror carefully with a lens tissue or lint-free cloth.

Sterilize the Scanner Tip

Scanner Tips shipped non-sterile. You must sterilize them before use.

CAUTION: For TP202 and TP203, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 60 cycles.

CAUTION: For TP204, if you limit the exposure time at 134°C to not more than 6 minutes, you can autoclave the Scanner Tip up to 180 cycles.

To sterilize the Scanner Tip, follow these steps:

- Put the Scanner Tip into a sealed steam sterilization pouch.
- Place the Scanner Tip in a steam autoclave for sterilization.
 - ◆ Exposure temperature should be set as 134°C.
 - ◆ Exposure time should exceed 3 minutes.
 - ◆ Exposure time should not exceed 6 minutes.

Precautions before Use

Perform the following activities on your product and accessories before use.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- ◆ Inspect the Lens Window.
- ◆ Inspect around the buttons and contact points.

If damage is noted, do not use the scanner and contact your representative.

Visually Inspecting the Scanner Tip

Visually inspect the Scanner Tip for signs of deterioration by doing the following:

- ◆ Verify that the Scanner Tip is not damaged, and its components are not detached.
- ◆ Verify that the Scanner Tip mirror does not have any smudges or scratches on it.

If deterioration is noted, replace the Scanner Tip.

WARNING

- The Lens Window on the scanner is a delicate optical component. Mount the lens window protective cover to protect the Lens Window from damage and dirt when the scanner is not in use.
- The mirror in the Scanner Tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.

In the event that you see poor scan quality or an unclear video preview in the software, clean the Scanner Tip mirror and the Lens Window of the Handpiece using a microfiber cleaning swab, applying ethanol that is free of impurities.

Troubleshooting

Problem Description	Action
The 3D model shows discrepancies and overlapping areas.	Perform the scan using the suggested scanning method. Use the Cut tool to eliminate mismatched data and excess tissue, then conduct a rescan.
Following the buccal bite registration, a gap or an overlap appear between the upper and lower jaws.	Adjust the Occlusal pressure level in Preferences. Delete the incorrect bite view and rescan.
A decline in precision or poorly stitched images is noticed during acquisition.	Clean the lens window at the scanner's base using a moist, lint-free cloth or lens tissue to remove dust or water stains. Ensure the Scanner Tip is securely attached, and check that the preview video does not have dark edges or other dirt.
Reconstructing metallic preparations can sometimes be challenging.	Adjust the scanner's position—such as its distance or angle—and capture more of the area during the scan. Move the surgical light away from the patient to decrease light scatter.
The Scanner Tip is installed but not detected. No live video is displayed, and the Scanner Tip not detected icon is displayed at the bottom-right of the interface.	After removing the scanner tip, wipe the Scanner Tip Detection Sensor to ensure it is free from dirt or debris. Reattach the Scanner Tip and verify that it is securely and firmly connected to the scanner.
The inner surface of the scanner's lens window at its base shows signs of fogging.	Attach a fully dry Scanner Tip to the scanner, place the scanner in the Holder or desk, and wait for the fogging to clear. If the fogging persists after 24 hours, contact your local service provider for help. Ensure the Scanner Tip is completely dry before attaching it, and avoid using disinfectant-soaked cloths to clean the scanner.
The scanner does not cast light, and a static preview image is displayed in the video preview screen.	Ensure that your computer meets the Minimum System Requirements. Close the software, and then launch the software again.
Scanner is not shown in the <i>Scanner List</i> in the software while the Status Indicator keeps blinking in breathing mode	Reset the network in the scanner. To do so, while the Status Indicator is in breathing mode, simultaneously press the Scan Button on both sides for two seconds and release. The Status Indicator should start blinking rapidly, then enter breathing mode. The device should then

	be shown in the <i>Scanner List</i> in the software.
--	--

Regulatory Information

Regulatory Information

The device complies with the following regulations:

MDR: (EU) 2017/745 Medical Device Regulation, Class I following the Rule 5.
FDA Center for Devices & Radiological Health CDRH - Title 21 CFR 872.3661 (USA).
Medical Devices Regulations (Canada).
RoHS: Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment, 2011/65/EU Annex II and its amendment Directive (EU) 2015/863
RED: Directive 2014/53/EU The Radio Equipment Directive
FCC: Part 15 of The Federal Communications Commission Rules
ISED: Innovation, Science and Economic Development Canada

Compliance with European and International Standards

EN / IEC 60601-1: Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance

ANSI/AAMI ES 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

CAN/CSA-C22.2 No. 60601-1: Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

EN / IEC 60601-1-2: Medical Electrical Equipment, Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests

EN / IEC 80601-2-60: Medical electrical equipment — Part 2-60: Particular requirements for the basic safety and essential performance of dental equipment

EN / IEC 62471: Photobiological safety of lamps and lamp systems

EN / ISO 17664: Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices

EN / ISO 17665-1: Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

EN / IEC 60601-1-6: Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability

EN / IEC 62366-1: Medical devices - Part 1: Application of usability engineering to medical devices

EN / IEC 62304: Medical device software - Software life cycle Processes

EN ISO 10993: Biological evaluation of medical devices

ISO 14971: Medical devices - Application of risk management to medical devices

EN / ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements

EN ISO 20417: Medical devices — Information to be supplied by the manufacturer

ISO 9687: Dentistry - Graphical symbols for dental equipment

AAMI TIR 12: Designing, testing and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers

AAMI TIR 30: A compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

EN / IEC 62133-2: Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems

EN 50566: Product standard to demonstrate the compliance of wireless communication devices with the basic restrictions and exposure limit values related to human exposure to electromagnetic fields in the frequency range from 30 MHz to 6 GHz: hand-held and body mounted devices in close proximity to the human body.

EN 301489-1: Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU and the essential requirements of article 6 of Directive 2014/30/EU

EN 301489-17: Electromagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems; Harmonized Standard covering the essential requirements of article 3.1(b) of Directive 2014/53/EU

EN 301893: 5GHz RLAN; Harmonized Standard covering the essential requirements of article 3.2 of Directive 2014/53/EU

Classification in Accordance with EN/IEC 60601-1

Type of protection against electric shock: Internally powered

Degree of protection against electric shock: Type BF Applied Part

Mode of operation: Continuous operation

Flammable anesthetics: Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/IEC 60601-1-2

IEC 60601-1-2 EMC requirements and tests, Medical Electrical Equipment including CISPR 11 Group 1, Class B.

Electromagnetic Compatibility Precautions

). Medical electrical devices need specific measures to ensure electromagnetic compatibility (EMC). Medical devices should be set up and used following the EMC guidelines detailed in this document.

Even equipment that meets CISPR emission standards may disrupt communication with the AS 200E.

WARNING: To maintain the AS 200 E's optimal performance, ensure that all portable and mobile RF communication devices—including peripherals like antenna cables and external antennas—are used at a distance of no less than 30 cm (12 inches) from any part of the system, including the cables specified by the manufacturer. Using such equipment closer than this recommended distance may result in decreased performance of the AS 200E.

WiFi

The device operates with the 802.11a/n/ac protocol. Channel 38 or 46 is priority used for the handpiece. The channel bandwidth is 40MHz. The frequency range is 5150-5250MHz, 5725-5850MHz (the actual frequencies are dependent on local regulations and the configuration of the product). The maximum output power is 17.88 dBm.

This device complies with part 15 of the FCC Rules and contains license exempt transmitter(s)/receiver(s) that comply with ISSED's license-exempt RSS(s).

Operation is subject to the following two conditions:

- This device may not cause interference.
- This device must accept any interference, including interference that may cause undesired operation of the device.

NOTE: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

SAR Value for Handpiece:

0.798 W/Kg, 10g for CE

0.011 W/Kg, 1g for FCC

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Limited by local law regulations, the version for North America does not have a region selection option.

The device is for indoor use only and operates in the 5150-5250 MHz band to reduce the potential for harmful interference to co-channel mobile satellite systems.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Restrictions in the 5 GHz band: According to Article 10 (10) of Directive 2014/53/EU, the packaging shows that this radio equipment will be subject to some restrictions when placed on the market in Belgium (BE), Bulgaria (BG), the Czech Republic (CZ), Denmark (DK), Germany (DE), Estonia (EE), Ireland (IE), Greece (EL), Spain (ES), France (FR), Croatia (HR), Italy (IT), Cyprus (CY), Latvia (LV), Lithuania (LT), Luxembourg (LU), Hungary (HU), Malta (MT), Netherlands (NL), Austria (AT), Poland (PL), Portugal (PT), Romania (RO), Slovenia (SI), Slovakia (SK), Finland (FI), Sweden (SE), Turkey (TR), Norway (NO), Switzerland (CH), Iceland (IS), and Liechtenstein (LI). The WLAN function for this device is restricted to indoor use only when operating in the 5150 to 5250 MHz frequency range.

AT, BE, BG, CH, CY, CZ, DE, DK, EE, EL, ES, FI, FR, HR, HU, IE, IS, IT, LI, LT, LU, LV, MT, NL, NO, PL, PT, RO, SE, SI, SK, TR

Guidance and Manufacturer's Declaration

Guidance and Manufacturer's Declaration - Electromagnetic Emission

The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1 Class B	The DEVICE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference

		in nearby electronic equipment.
--	--	---------------------------------

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems (from Non-RF wireless communication equipment)

The AS 200E should be operated only within the electromagnetic environments described below (IEC 60601-1-2) for optimal performance. Customers or users must ensure that the device is used under these specified conditions.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	For optimal conditions, floors should be constructed from wood, concrete, or ceramic tile. If synthetic materials are used to cover the floors, please ensure that the relative humidity is maintained at a minimum of 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For optimal conditions, power-frequency magnetic fields should be maintained at levels commonly found in standard commercial or hospital environments.
Radiated RF IEC 61000-4-3	3 V/m 80MHz – 2.7GHz	3 V/m 80MHz – 2.7GHz	Environment of a professional healthcare facility.

NOTE: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEVICE is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the DEVICE.

Guidance and Manufacturer's Declaration – Electromagnetic Immunity for Equipment and Systems (from RF Wireless Communication Equipment)

The AS 200E meets the immunity requirements for proximity fields from RF wireless communication devices, in line with the test levels outlined below according to the IEC 60601-1-2 standard. It is the customer or user's responsibility to ensure that the AS 200E is operated in a suitable environment.

Test Frequency (MHz)	Band (MHz)	Immunity Test Levels
385	380 – 390	Pulse modulation 18Hz, 27V/m
450	430 – 470	FM, ± 5 kHz deviation, 1 kHz sine, 28V/m
710	704 – 787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800 – 960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700 – 1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400 – 2570	Pulse modulation 217Hz, 28V/m
5240	5100 – 5800	Pulse modulation 217Hz, 9V/m
5500		
5785		

Accessories

Using cables or accessories other than those specified—unless they are manufacturer-approved replacements for internal components—could lead to increased emissions or reduced immunity of the medical device.

Other Equipment

WARNING: Avoid placing this equipment next to or stacking it with other devices, as this may cause improper functioning. If such placement is necessary, monitor this equipment and the other devices to ensure they operate correctly.

Technical Specifications

Product Model

AS 200E

Technical Specifications

Components	Technical Specifications
Weight	Scanner (including battery, without Scanner Tip): 245 g Tip Large - Sterilization (TP202): 14g Tip Small - Sterilization (TP203): 12g Tip Large - Sterilization (TP204): 14g
Color	3D Full Color
USB Wireless Card Port	USB 3.0
USB Charge Port	Type C Port
Field of View	Tip Large - Sterilization (TP202): 16mm x 14mm Tip Small - Sterilization (TP203): 12mm x 12mm Tip Large - Sterilization (TP204): 16mm x 14mm
Depth of Field	15mm
Minimum System Requirements of Laptop	Processor: 10th Generation Intel® Core™ i7-10750H or AMD Ryzen™ 7 4700U Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen

Minimum System Requirements of Desktop	Processor: 9th Generation Intel® Core™ i5-9600 or AMD Ryzen™ 5 2600 Memory: 16 GB Graphics Card: NVIDIA GeForce GTX 1650 4GB Disk: 512 GB SSD Operating System: Win 10 (build 18362+) / Win 11, 64 bit Others: USB 3.0 Port Optional: Touch Screen
---	---

CAUTION: It is MANDATORY to check that your system configuration is compatible with the computer system requirements for the accompanying software.

Environmental Requirements

Components	Environmental Requirements
Operating Temperature	15°C ~ 30°C
Transport and Storage Temperature	-10°C ~ 50°C
Operating Relative Humidity	10% ~ 65% RH
Transportation and Storage Relative Humidity	10% ~ 95% RH
Operating Atmospheric Pressure	70 ~ 106 KPa
Transportation and Storage Atmospheric Pressure	60 ~ 106 KPa

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Etkon GmbH Lochhamer Schlag 6, 82166 Gräfelfing, Germany