
ScanPro

Instructions, Vejledning,
Gebrauchsanweisung, Instrucciones,
Instructions, Istruzioni, Instructies,
Bruksanvisning, Instruktioner, 使用説明書

Contents

4 EN/English

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EN/English

» Notice

The information in this document is subject to change. Neither Alliedstar nor any of its subsidiaries shall be liable for errors contained herein or for incidental damages in conjunction with the furnishing, performance, or use of this material. No part of this publication may be reproduced without the permission of Alliedstar.

We recommend that you thoroughly familiarize yourself with this Guide to make the most effective use of your system.

All trademarks and registered trademarks are the property of their respective holders.

U.S. Federal law restricts this device to sale by or on the order of a dentist.

Conventions

The following special messages emphasize information or indicate potential risks to personnel or equipment.


	WARNING: Avoid injury to yourself or others by following the safety instructions precisely.
	CAUTION: Alerts you to a condition that might cause serious damage or cause problems.
	NOTE: Provides extra information and hints.

» Software Overview

Interface

The scanner operates with ScanPro software. The software interface enables you to acquire 3D models in two ways:

- Partial arch scan: Several teeth in the preparation area on both the upper and lower jaws, and buccal bite registration.
- Full arch scan: Upper jaw, lower jaw, and buccal bite registration.



01	Acquisition Step: Displays the current step in the acquisition process.
02	Jaw/Bite Switch: Enables you to select the upper jaw, lower jaw, or buccal bite registration.
03	Image Toolbar: Enables you to select and manage 3D model.
04	Soft tissue removal Switch: Enables you to automatically remove soft tissue.
05	Shining surface: When scanning the highlight metal, you can turn on the button to scan.
06	Impression Scan: Turn on this function to scan an impression, after the scanning is completed, the model will be inverted and merged into existing model.
07	Workflow Toolbar: Enables you to choose general scan, preparation scan and scan body scan.
08	3D Model Display Screen: Displays the 3D model created from the scan.
09	Video Preview Screen: Displays live video when scanning, or the scanner status when not scanning.
10	Scanning Time: Displays the total scanning duration of the current case.
11	Display Toolbar: Enables you to change how the 3D model is displayed.
12	Option Menu: Access the system configuration and version information.
13	AS Connect Service: You can view the cloud service status.
14	Login Menu: Enables you to log in, switch users and log out.

15	Battery: Displays the battery level of wireless scanner, if a wired scanner is connected this icon will not present.
----	---

Toolbar

Acquisition Step Toolbar	
 Scan	Scan Button: Enables you to scan the upper and lower arch, and the buccal bite registration.
 Refine	Refine button: Refines the acquired 3D model and enables you to use various tools to check the refined results.
 Finish	Finish button: complete the case information and save the scan results.
Jaw/Bite Switch Toolbar	
	Upper Jaw button: Acquires a 3D model of the upper jaw.
	Lower Jaw button: Acquires a 3D model of the lower jaw.
	Buccal Bite Registration button: Acquires a 3D model of the bite.
	Switch Upper/Lower Jaw button: Changes the acquisition mode from upper to lower or vice versa, if you accidentally scan teeth on the wrong jaw.
Image Toolbar	
	Cut button: draw a curve to delete unnecessary data.
	Undo last cut
	Return to the upper-level (icon on the left)
	Lock button: Enables you to select and lock an area on the model to prevent it from being updated by additional scanning.
	Unlock the last locked area
	Return to the upper level (icon on the left)
	Mark tooth button: Mark one or more preparation / implant areas. You can use the wheel on your mouse to adjust the diameter of the mark.
	Delete marked teeth
	Return to the upper level (icon on the left)
	Delete button: Delete all models from the current case.

	Add impression region: Choose one or more regions on the model, and only align data in the selected regions.
	Unlock the last locked area
	Return to the upper level (icon on the left)
	Quadrant Snapshot button: Displays a preview of five 2D images showing different views of the model.
 	Transparency button: Set the transparency of the scanned model.
	Undercut check button: View the undercut areas on the surface of the preparation.
	Avoid Intersections: When this function is turned on, software will automatically adjust the bite registration to avoid bite intersections. This function is turned off by default.
	Occlusion analysis button: Analyze the occlusion space.
	Measurement button: Measure the distance from the specified area of the section.
	Margin line button: Automatically or manually create margin lines, and edit margin lines.
Display Toolbar	
	True Color button: When selected, displays the 3D model in the actual color of the patient's hard and soft tissues. When deselected, displays the 3D model in monochrome.
	Scan Quality Detection button: Turn on or turn off color overlay that indicates the quality of the scan. Regions with better scan quality will be green.
	Show 3D center button: When selected, the rotation center will be displayed when rotating the 3D model.
	Zoom fit button: Scales the 3D model to the best size to fit the display region
	View orientation button: Displays a list of views.
	Front view
	Back view

	Left view
	Right view
	Top view
	Bottom view
	Return to the up-level toolbar (icon on the right)
	Snapshot button: Takes a snapshot of the 3D model as it appears on your screen.
	Intraoral Camera button: Enables you to select 2D intraoral images
Workflow Toolbar	
	Common Scan Button: Enables you to acquire a 3D model of the upper jaw, lower jaw, and the buccal bite.
	Config scan button: Edit acquisition type.
	Add preparation scan button: Add a preparation scan step.
	Add scan body button: Add a scan body scan step.
	Preparation button: Activate the preparation scan.
	Scan body button: Activate the scan body scan.
Option Menu	
	Option menu button: Enables you to access Import/Export scan data, Scan history, Preference and About dialog.
Scanner Status Icons	
	Scanner is connecting
	Scanner Tip is not detected
	Scanner is overheating

	Scanner is in pause state
	Scanner is in sleep mode
	Scanner is in the (charging) holder

Option Menu

Import/Export Scan Data

The import/export scan data menus allow users to export the current scan data to a specified folder, and import it later to continue scanning or perform other operations.

CAUTION: If you need to continue scanning after importing the previously saved scan data, you must ensure that the scan data is acquired by the same scanner currently connected, otherwise you will not be able to perform subsequent scans on the imported data.

Scan History Dialog

The Scan History dialog allows users to import automatically saved scan history data. Based on different user settings, the scan history data is retained for up to 30 days. In the scan history dialog, you can search data by ID or name. To import or manage the selected record, right-click the entry, and then select “Open” , “Lock/Unlock” or “Delete” .

Open: Import the selected scan history data.

Lock/Unlock: Lock or unlock the selected scan history. The locked scan history records will not be cleared after the storage period expires.


Delete: Delete the selected scan history data.

CAUTION: If you need to continue scanning after importing the previously saved scan data, you must ensure that the scan data is acquired by the same scanner currently connected, otherwise you will not be able to perform subsequent scans on the imported data.

Save To Scan History: If users need to save scan data to Scan history after scanning, you can click this button.

Preference Dialog

The Preferences dialog enables you to select the software and scanner settings.


Language: Select the user interface language.

Tooth numbering system: Select the FDI/ISO or Universal as the tooth numbering system.

Automatically update software and firmware: When selected, the software will automatically prompt of new versions.

Show patient information dialog on startup: When selected, the patient information dialog will pop up when the software opens.

Show instructions: When selected, the software will show tutorial animation.


Hole highlight color: Fill the holes in the model with the specified color after refining.

Enable remote control: When selected, the scanned 3D model can be rotated by rotating the scanner after the scanning is stopped.

Enable scanning sound: When selected, a sound will be continuously played when you are successfully scanning (if your computer does not have speakers, this option will not take effect).

Enable warning sound: When selected, a warning sound will be played if the scanning duration exceeds the recommended thresholds, a strong light is detected, or the scanning performance is declined (if your computer does not have speakers, this option will not take effect).


Color scheme: Select the color scheme of the user interface.

UI element size: Select the size of interface icons to adapt to different screen resolutions.

Monochrome color: Select the color when displaying the 3D model in monochrome.

Jaw/Bite switch: Select the style of jaw/bite switch.


Save bite data: When selected, save the bite images to separate files.

HD orthodontic data: When this option is selected, scanning data will be saved in higher fidelity.

HD color: When this option is selected, texture will be saved in a separate file with higher fidelity.

Automatically save scan history: When selected, scan history is automatically saved when the software is closed. When this option is enabled, the user can customize the number of days and the path to save the scan history.


Shut down the wired scanner after: User can customize the automatic shutdown time of the wired scanner.

Shut down the wireless scanner after: User can customize the automatic shutdown time of the wireless scanner.

Create support file: Export software and firmware logs.

About Dialog

The About dialog displays software version, firmware version and scanner serial number information.

» Acquiring a 3D Model Using the Common Scan Workflow

The common scan workflow (default workflow) allows you to scan either a partial or full arch. For either type of scan, you should scan upper jaw, lower jaw, and buccal bite registration. Once this basic information has been scanned, other scan types (such as scan body scanning) can be used to acquire other information.

For example, use a general scanning workflow to scan a 3D model for orthodontics. Scan the complete dental arch to create a 3D model. The user should scan the upper jaw, lower jaw, and buccal bite registration, as well as images of a few millimeters of gum tissue in the lingual and buccal scans. The software will combine these images to create a 3D model, which can then be uploaded to AS Connect or saved locally.

NOTE: In some cases, you can acquire 3D models of a single arch (partial or full) and not obtain a buccal bite registration (for example, if there are no teeth in the opposing arch), but it is recommended that you acquire both arches and a buccal bite registration when possible.


To acquire a 3D model, follow these steps:

- Scanning the upper and lower jaw.
- Scan the buccal bite registration.
- Refine and check the 3D model.
- Complete and save the 3D model.


Scanning the Upper and Lower Jaw

To scan a 3D model of the upper and lower jaw, follow these steps:

1. Dry the teeth thoroughly before starting an acquisition.
2. In the interface, select the Upper Jaw acquisition mode  or press the mode button on the scanner to select the upper jaw scan mode.
3. Place the Scanner Tip of the scanner on the surface of the tooth to stabilize the scanner and press the Scan Button. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-5mm from the teeth. The image will be automatically scanned and displayed on the 3D model display area.

4. Slowly move the Scanner Tip of the scanner along the occlusal surface to scan the remaining teeth in the arch.

NOTE: The holes on the 3D model will be displayed in the color specified by the user. It is recommended to scan these areas until the holes disappear.


CAUTION: Re-dry the teeth as appropriate throughout the acquisition process.

5. When the occlusal surface scan is complete, scan the lingual or buccal surface of the teeth in the arch.

6. When the lingual or buccal scan is completed, scan the opposite side of the arch.

If you need to remove soft-tissue artifacts, mismatches, or unwanted views during the acquisition, click the Cut tool , then draw a curve that covers the region to be removed on the 3D model. If necessary, rescan the area to fill the holes.

NOTE: You can click and hold the middle mouse button to rotate the 3D model when using tools such as Cut , Lock  and Mark tooth .

7. After scanning the upper jaw, you can continue scanning the lower jaw. In the interface, select the Lower Jaw acquisition mode  or press the mode button on the scanner to select the lower jaw scan mode.

8. Repeat steps 3 to 6 until the lower jaw scan is completed.

9. Check whether there are obvious holes on the 3D model and rescan if necessary.

10. When you confirm that the 3D model of the upper and lower jaw is complete, continue with the buccal bite registration acquisition.

Scanning the Buccal Bite Registration

To acquire a buccal bite registration, follow these steps:

- In the interface, select the **Buccal Bite Registration**  acquisition mode or press the mode button on the scanner to select the buccal bite registration mode.
- Place the Scanner Tip of the scanner into the buccal side in the patient's mouth, then rotate the Scanner Tip to align with the teeth, close the patient's mouth and confirm that the bite position is correct.
- Press the Scan Button on the scanner, slowly move the Scanner Tip in mesial direction with equal coverage of upper and lower teeth. The example below shows a buccal bite registration. You can use the toolbar on the left side to switch between the scanned occlusion views or delete the scanned occlusion views for rescanning.
- Turn on the **Avoid intersections** function to automatically adjust the bite registration to avoid bite intersections.

- Use Occlusion analysis to analyze the occlusion space.


NOTE: You can scan one or two buccal bite registrations. It is recommended to scan one on the left side and one on the right side of the patient's mouth.

- After scanning the buccal bite registrations, rotate the model and zoom the view to ensure that the bite is accurate and that there are no areas where the bite is mismatched. If necessary, you can delete the scanned occlusion and rescan.

Refining and Checking the 3D Model

Refining the 3D model allows you to obtain higher accuracy data for further processing. To refine the 3D model, follow these steps:

- Click the **Refine** button, the refining progress bar will be displayed. Depending on your computer configuration, the refining process may take several minutes.


- After the refining is complete, manipulate the refined 3D model using the following methods:
 - » Click and hold on the 3D model to rotate it.
 - » Right-click and hold on the 3D model to move it in the window.
 - » Use the wheel of the mouse to zoom in or zoom out on the 3D model.
 - » Click the **True color** button to view the color or monochrome 3D model.

- » Click the **Zoom fit** button  to scale the model to its best view.
- » Click the **View orientation** button  to see six views of the 3D model.
- » Click the **Snapshot** button  to take a snapshot of the 3D model.
- » Click **Upper Jaw** button  or **Lower Jaw** button  to hide or display the jaws.
- » Click the **Transparency** button  to display the transparency slider, click and drag the slider to adjust the transparency of the 3D model.
- » Click **Cut** button  to select and delete unwanted data.
- » Click Intraoral **Camera** button  to pick up intraoral images from the scanned data.
- » Click the **Quadrant Snapshot** button  to open the Quadrant Snapshot window, which displays multiple views of the model.
- » Turn on the **Avoid intersections**  function to automatically adjust the bite registration to avoid bite intersections.
- » Click the **Occlusion analysis** button  to analyze the occlusion space.
- » Click the **Measurement** button  to measure the distance between two points on a clipping plane.
- If you find obvious holes when checking the 3D model, click the **Scan Button**  **Scan** and rescan the areas to fill the holes.
- Repeat these steps until you are satisfied with the 3D model.

Completing and Saving the 3D Model

To complete the scan and save the 3D model, follow these steps:

- Click the **Finish** button  **Finish** and the following page will be displayed. If you entered patient information in the previous step, it will be displayed on this page.

- Complete the patient information, if necessary, you can add some attachments to the case.
- Click the **Save to file system** button  to save the case and 3D model data to the specified local directory.

» Acquiring a 3D Model Using the Scan Body Workflow

You can use the scanner to scan all or part of the dental arch containing the abutment or implant and create a 3D model. When using an existing abutment to scan all or part of the dental arch, use the “common scan” workflow. When you scan all or part of the dental arch that contains the implant, the scan body workflow is used. When using the scan body workflow, the software will copy the 3D model and automatically cut out the implant areas, rescan the tooth jaw with the scan body and create two 3D models: one contains the scan body and the other does not contain the scan body. The user can send these models to the dental laboratory.

NOTE: In some cases, you can acquire 3D models of a single arch (partial or full) and not obtain a buccal bite registration (for example, if there are no teeth in the opposing arch), but it is recommended that you acquire both arches and a buccal bite registration when possible.


To acquire a 3D model for an implant, follow these steps:

- Scan the upper and lower jaw.
- Scan the buccal bite registration.
- Mark the implant areas.
- Install and scan the scan body.
- Refine and check the 3D model.
- Complete and save the 3D model.

Scanning the Upper and Lower Jaw

To scan a 3D model of the upper and lower jaw, follow these steps:

11. Dry the teeth thoroughly before starting an acquisition.
12. In the interface, select the Upper Jaw acquisition mode  or Press the mode button on the scanner to select the upper jaw scan mode.
13. Place the Scanner Tip of the scanner on the surface of the tooth to stabilize the scanner and press the **Scan Button**. Wait until a 3D model appears on the 3D model display screen, and then slowly move it along the arch at 0-5mm from the teeth. The image will be automatically scanned and displayed on the 3D model display area.

14. Slowly move the Scanner Tip of the scanner along the occlusal surface to scan the remaining teeth in the arch.


NOTE: The holes on the 3D model will be displayed in the color specified by the user. It is recommended to scan these areas until the holes disappear.

15. When the occlusal surface scan is complete, scan the lingual or buccal surface of the teeth in the arch.

16. When the lingual or buccal scan is completed, scan the opposite side of the arch.

CAUTION: Re-dry the teeth as appropriate throughout the acquisition process.

If you need to remove soft-tissue artifacts, mismatches, or unwanted views during the acquisition, click the Cut tool , then draw a curve that covers the region to be removed on the 3D model. If necessary, rescan the area to fill the holes.

17. After scanning the upper jaw, you can continue scanning the lower jaw. In the interface, select the Lower Jaw acquisition mode  or press the mode button on the scanner to select the lower jaw scan mode.

18. Repeat steps 3 to 6 until the lower jaw scan is completed.

19. Check whether there are obvious holes on the 3D model and rescan if necessary.
20. When you confirm that the 3D model of the upper and lower jaw is complete, continue with the buccal bite registration acquisition.

Scanning the Buccal Bite Registration

To acquire a buccal bite registration, follow these steps:

- In the interface, select the **Buccal Bite Registration**  acquisition mode or press the mode button on the scanner to select the buccal bite registration mode.
- Place the Scanner Tip of the scanner into the buccal side in the patient's mouth, then rotate the Scanner Tip to align with the teeth, close the patient's mouth and confirm that the bite position is correct.
- Press the Scan Button on the scanner, slowly move the Scanner Tip in mesial direction with equal coverage of upper and lower teeth. The example below shows a buccal bite registration. You can use the toolbar on the left side to switch between the scanned occlusion views, or delete the scanned occlusion views for rescanning.
- Turn on the Avoid intersections function to automatically adjust the bite registration to avoid bite intersections.
- Use Occlusion analysis to analyze the occlusion space.

NOTE: You can scan one or two buccal bite registrations. It is recommended to scan one on the left side and one on the right side of the patient's mouth.

- After scanning the buccal bite registrations, rotate the model and zoom the view to ensure that the bite is accurate and that there are no areas where the bite is mismatched. If necessary, you can delete the scanned occlusion and rescan.

Marking the Implant Areas

After the scan is complete, you can mark the implant areas so that you can rescan the area after installing the scan body. To mark one or more implant areas, follow these steps:

- In the interface, select the jaw  with implant to activate the 3D model.
- Click the Mark tooth button , and then click in the center of the implant. A 3D ball will display on the occlusal surface indicating an implant. To adjust the diameter of the cylinder, hover the cursor over it and scroll the wheel.

- Rotate the 3D model if necessary, and mark all the implant areas.

NOTE: You can click and hold the middle mouse button to rotate the 3D model when using tools such as Cut , Lock  and Mark tooth .

- Click the Return button () to quit Mark tooth tool.

Installing and Scanning the Scan Body

After installing the scan body, follow the steps below to scan the scan body:


- After making the implant areas, click the **Config Scan Button** at the bottom of the interface.


- Click the **Add scan body** button in the pop-up window.

- Click the **Scan Body** button below.

- Check whether the part that the software automatically cut out from the implant areas is sufficient to contain the diameter of the scan body. If the cut out is not enough to contain the scan body, first click Common Scan icon, then click Mark tooth button in the Mark Tooth section and use mouse wheel to adjust the diameter of the cut out section. Then click the Config Scan button and select Add scan body to apply the change.
- Start the scanner and scan the part of the scan body areas that has been cut out.


NOTE: Start the scan body scanning from the surface 1-2 teeth away from the scan body, so that the system can recognize the 3D structure.

- When you have completed the scan body scanning, proceed to the **Refining and Checking the 3D Model** step.

Refining and Checking the 3D Model

Refining the 3D model allows you to obtain higher accuracy data for further processing. To refine the 3D model, follow these steps:

- Click the **Refine** button  **Refine** button, the refining progress bar will be displayed. Depending on your computer configuration, the refining process may take several minutes.


- After the refining is complete, manipulate the refined 3D model using the following methods:
 - » Click and hold on the 3D model to rotate it.
 - » Right-click and hold on the 3D model to move it in the window.
 - » Use the wheel of the mouse to zoom in or zoom out on the 3D model.
 - » Click the **True color** button  to view the color or monochrome 3D model.
 - » Click the **Zoom fit** button  to scale the model to its best view.

- » Click the **View orientation** button  to see six views of the 3D model.
- » Click the **Snapshot** button  to take a snapshot of the 3D model.
- » Click **Upper Jaw** button  or **Lower Jaw** button  to hide or display the jaws.
- » Click the **Transparency** button  to display the transparency slider, click and drag the slider to adjust the transparency of the 3D model.
- » Click **Cut** button  to select and delete unwanted data.
- » Click Intraoral **Camera** button  to pick up intraoral images from the scanned data.
- » Click the **Quadrant Snapshot** button  to open the Quadrant Snapshot window, which displays multiple views of the model.
- » Turn on the Avoid intersections  function to automatically adjust the bite registration to avoid bite intersections.
- » Click the **Occlusion analysis** button  to analyze the occlusion space.
- » Click the **Measurement** button  to measure the distance between two points on a clipping plane.
- If you find obvious holes when checking the 3D model, click the **Scan Button**  **Scan** and rescan the areas to fill the holes.
- Repeat these steps until you are satisfied with the 3D model.

Completing and Saving the 3D Model

To complete the scan and save the 3D model, follow these steps:

- Click the **Finish** button  **Finish** and the following page will be displayed. If you entered patient information in the previous step, it will be displayed on this page.

- Complete the patient information, if necessary, you can add some attachments to the case.
- Click the **Save to file system** button  to save the case and 3D model data to the specified local directory.

» Acquiring a 3D Model Using the Preparation Workflow

Depending on your workflow preferences, you can use the preparation workflow. This workflow provides you with two choices:

- You can scan the tooth before preparing, and then take another scan after the tooth is prepared.
- You can import the previously saved data before preparing, and then start to scan the preparation.


With either option, two 3D models will be created after the scan is completed: one containing the unprepared tooth and the other containing the preparation. You or the lab can use these 3D models for restorations to make them closely resembles the original tooth.

The following steps describe how to combine the previous scan with the new scan of the preparation.

- Import the scanned data before preparing.
- Scan the preparation areas.
- Refine and check the 3D model.
- Complete and save the 3D model.

Import the scanned data before preparing

21. Click the **Option menu**  button and select Import scan data.
22. Select the .scan file you want to import, and click Open. Wait for the software to import the data, and the 3D model will be displayed after the import is complete.

- Click the Mark tooth button , then click in the center of the occlusal surface of the preparation. A 3D ball will display on the occlusal surface indicating an implant. To adjust the diameter of the cylinder, hover the cursor over it and scroll the wheel.
- Rotate the 3D model if necessary, and mark all the preparation areas.


NOTE: You can click and hold the middle mouse button to rotate the 3D model when using tools such as Cut , Lock  and Mark tooth .

- Click the Return button () to quit Mark tooth tool


Scan the Preparation Areas

After marking the preparation areas, you need to rescan these areas. To scan one or more marked preparation areas, follow these steps:


- Click the **Config Scan Button** at the bottom of the interface.


- Click the **Add preparation button** in the pop-up window.

- Click the **Preparation** button below

- Check whether the part that the software automatically cut out from the implant areas is sufficient to contain the diameter of the scan body. If the cut out is not enough to contain the scan body, first click Common Scan icon, then click Mark tooth button in the Mark Tooth section and use mouse wheel to adjust the diameter of the cut out section. Then click the Config Scan button and select Preparation to apply the change.
- Start the scanner and scan the part of the preparation areas that has been cut out.

- If you need to remove soft-tissue artifacts, mismatches, or unwanted views during the acquisition, click the **Cut** tool , then draw a curve that covers the region to be removed on the 3D model. If necessary, rescan the area to fill the holes.

NOTE: You can click and hold the middle mouse button to rotate the 3D model when using tools such as Cut , Lock  and Mark tooth .


- When you have completed the preparation scanning, proceed to the **Refining and Checking the 3D Model**

step.

Refining and Checking the 3D Model

Refining the 3D model allows you to obtain higher accuracy data for further processing. To refine the 3D model, follow these steps:

- Click the **Refine**  **Refine** button, the refining progress bar will be displayed. Depending on your computer configuration, the refining process may take several minutes.


- After the refining is complete, manipulate the refined 3D model using the following methods:
 - » Click and hold on the 3D model to rotate it.
 - » Right-click and hold on the 3D model to move it in the window.
 - » Use the wheel of the mouse to zoom in or zoom out on the 3D model.
 - » Click the **True color** button  to view the color or monochrome 3D model.
 - » Click the **Zoom fit** button  to scale the model to its best view.
 - » Click the **View orientation** button  to see six views of the 3D model.
 - » Click the **Snapshot** button  to take a snapshot of the 3D model.
 - » Click **Upper Jaw** button  or **Lower Jaw** button  to hide or display the jaws.
 - » Click the **Transparency** button  to display the transparency slider, click and drag the slider to adjust the transparency of the 3D model.
 - » Click **Cut** button  to select and delete unwanted data.
 - » Click Intraoral **Camera** button  to pick up intraoral images from the scanned data.
 - » Click the **Quadrant Snapshot** button  to open the Quadrant Snapshot window, which displays multiple views of the model.
 - » Turn on the Avoid intersections  function to automatically adjust the bite registration to avoid bite intersections.
 - » Click the **Occlusion analysis** button  to analyze the occlusion space.
 - » Click the **Measurement** button  to measure the distance between two points on a clipping plane.
- If you find obvious holes when checking the 3D model, click the **Scan Button**  **Scan** and rescan the areas to fill the holes.
- Repeat these steps until you are satisfied with the 3D model.

Completing and Saving the 3D Model

To complete the scan and save the 3D model, follow these steps:

- Click the **Finish** button  and the following page will be displayed. If you entered patient information in the previous step, it will be displayed on this page.

The screenshot displays the ScanPro application interface during the 'Finish' step of the 3D model preparation workflow. The top navigation bar includes the 'ScanPro' logo, a 99% battery indicator, and a user profile icon. The main workflow progress bar shows 'Scan', 'Refine', and 'Finish' (the active step). The central panel is divided into three sections: 'Patient Information' with fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment box; 'Tooth Selection' showing a dental arch diagram with a blue figure; and 'Attachments' with a plus sign icon. A 'Voice Memo' section at the bottom right includes a microphone icon, a play button, and a delete icon. On the far right, there are icons for uploading and downloading files.

- Complete the patient information, if necessary, you can add some attachments to the case.
- Click the **Save to file system** button  to save the case and 3D model data to the specified local directory.

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
The screenshot shows the ScanPro application interface. At the top, there's a status bar with 'ScanPro' and a battery icon at 99%. Below it, a progress bar shows three steps: 'Scan' (disabled), 'Refine' (disabled), and 'Finish' (active, highlighted in blue). The main content area is divided into three panels:

- Patient Information:** Contains fields for ID (001), Name (John Doe), Gender (Male/Female radio buttons), Age, and a Comment box (0/200 characters).
- Tooth Selection:** Displays a circular diagram of teeth with a blue figure of a person inside, indicating the patient's position.
- Attachments:** A large area with a plus sign icon for adding attachments.

On the right side of the main content area, there are two icons: a cloud upload icon and a download icon. Below these is a 'Voice Memo' section with a microphone icon, a play button, and a delete icon.

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
The screenshot shows the 'Finish' step of the ScanPro software. The top bar has three buttons: 'Scan', 'Refine', and 'Finish'. The main area is divided into three panels:

- Patient Information:** Contains fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and Comment (0/200).
- Tooth Selection:** Shows a diagram of teeth with a blue figure indicating a selection.
- Attachments:** Includes a plus sign for adding attachments and a 'Voice Memo' section with a microphone icon.

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The screenshot shows the ScanPro software interface during the 'Finish' step. The top navigation bar includes 'Scan', 'Refine', and 'Finish' buttons. The main content area is divided into three panels:

- Patient Information:** Contains fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment box (0/200).
- Tooth Selection:** Displays a circular diagram of teeth with a blue figure in the center.
- Attachments:** Features a plus sign icon for adding attachments.
- Voice Memo:** Includes a microphone icon, a play button, and a delete icon.

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
The screenshot shows the 'Finish' step of the ScanPro software. The top navigation bar includes 'Scan', 'Refine', and 'Finish' buttons. The main window displays a 'Patient Information' form with the following fields: ID (001), Name (John Doe), Gender (Male/Female), Age, and Comment. To the right of the form is an 'Attachments' section with a plus icon. At the bottom right is a 'Voice Memo' section with a microphone icon. The status bar at the top shows 'ScanPro' and '99%'.

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
The screenshot shows the ScanPro application window. At the top, there's a status bar with 'ScanPro' and a 99% battery icon. The main window has a header with 'User' and system icons. Below the header is a navigation bar with three buttons: 'Scan' (disabled), 'Refine' (disabled), and 'Finish' (active). The 'Finish' button is highlighted with a checkmark icon. The main content area is divided into three panels. The left panel, titled 'Patient Information', contains fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment box. The middle panel, titled 'Tooth Selection', shows a circular diagram of teeth with a blue figure of a person. The right panel, titled 'Attachments', has a large plus sign icon. At the bottom right, there's a 'Voice Memo' section with a microphone icon, a play button, and a delete icon.

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
- After the refining is complete, manipulate the refined 3D model using the following methods:
 - » Click and hold on the 3D model to rotate it.
 - » Right-click and hold on the 3D model to move it in the window.
 - » Use the wheel of the mouse to zoom in or zoom out on the 3D model.
 - » Click the **True color** button to view the color or monochrome 3D model.
 - » Click the **Zoom fit** button to scale the model to its best view.

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The screenshot shows the ScanPro software interface. At the top, there is a status bar with 'ScanPro' and a battery icon at 99%. Below this is a navigation bar with three buttons: 'Scan', 'Refine', and 'Finish'. The main area displays a form for 'Patient Information' and 'Attachments'. The 'Patient Information' section includes fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment field (0/200). The 'Attachments' section has a large plus sign icon for adding files. To the right of the form, there are two icons: a cloud upload icon and a download icon. At the bottom right, there is a 'Voice Memo' section with a microphone icon, a play button, and a delete button.

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The screenshot shows the ScanPro software interface during the 'Refine' step. The top bar includes the 'ScanPro' logo, a battery level indicator at 99%, and a user profile icon. The main navigation bar has three buttons: 'Scan', 'Refine' (active), and 'Finish'. The central workspace is divided into three panels:

- Patient Information:** Contains input fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment field (0/200 characters).
- Tooth Selection:** Displays a circular diagram of teeth with a blue figure indicating a selection.
- Attachments:** Features a large plus sign for adding attachments and a 'Voice Memo' section with a microphone icon, a play button, and a delete icon.

On the right side of the interface, there are two icons: an upload icon (cloud with an arrow) and a download icon (cloud with a downward arrow).

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
The screenshot displays the ScanPro application interface during the 'Finish' step of a 3D model preparation workflow. The top navigation bar includes 'Scan', 'Refine', and 'Finish' (highlighted). The main content area is divided into three panels:

- Patient Information:** Contains input fields for ID (001), Name (John Doe), Gender (Male/Female), Age, and a Comment field (0/200).
- Tooth Selection:** Features a circular diagram of teeth with a blue icon of a person and a dog in the center.
- Attachments:** Includes a large plus sign icon for adding attachments and a Voice Memo section with a microphone icon, a play button, and a delete icon.

On the right side of the interface, there are two icons: a cloud with an up arrow and a cloud with a down arrow.

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Copyright © 2022 Alliedstar
All rights reserved